

## Plural Word Lists for the Validation of Frequency Lowering

### Instructions for live voice validation

With frequency lowering deactivated, complete the following steps:

1. With mouth obscured, say a single word at a conversational level.
2. Ask the patient to repeat the word.
3. Mark the response as correct or incorrect, completing a single list.
4. Calculate the patient's score (% correct).

With frequency lowering activated, use an alternate list and repeat steps 1-4. If the patient's ability to detect the plural cues does not improve, increase Spectral iQ Bandwidth and Gain. Continue this process until the patient achieves 50 percent correct plural detection.

LIST 1	SCORE	LIST 2	SCORE	LIST 3	SCORE
Key		Play		Coats	
Days		Ears		Floor	
Plays		Keys		Signs	
Toy		Clouds		Eye	
Cities		Brick		Book	
Bricks		City		Frogs	
Pig		Way		Boats	
Cloud		Toys		Card	
Ways		Pigs		Caves	
Ear		Day		Pen	
TOTAL		TOTAL		TOTAL	

